


NanoJ: High-Performance Open-Source Super-Resolution Microscopy in ImageJ


R.F. Laine^{1,2*}, K.L. Tosheva^{1*}, R.D.M. Gray¹, P. Almada¹, D. Albrecht¹, J. Mercer¹, C. Leterrier³, P.M. Pereira^{1,2}✉, S. Culley^{1,2}✉ and R. Henriques^{1,2}✉; *equal contribution

¹University College London, UK; ²The Francis Crick Institute, UK; ³Aix-Marseille Univ, CNRS, INP, France
Correspondence to: p.pereira@ucl.ac.uk, s.culley@ucl.ac.uk, r.henriques@ucl.ac.uk

Poster DOI: 10.6084/m9.figshare.7963556
github.com/HenriquesLab — @HenriquesLab




Our team has built an open-source image analysis framework for super-resolution microscopy designed to combine high performance and ease of use. We named it NanoJ - a reference to the popular ImageJ software it was developed for. Here we highlight the current capabilities of NanoJ for several essential processing steps: spatio-temporal alignment of raw data (NanoJ-Core), super-resolution image reconstruction (NanoJ-SRRF), image quality assessment (NanoJ-SQUIRREL), structural modelling (NanoJ-VirusMapper) and control of the sample environment (NanoJ-Fluidics). We expect to expand NanoJ in the future through the development of new tools designed to improve quantitative data analysis and measure the reliability of fluorescent microscopy studies.



NanoJ-Core: Drift Correction


Drift estimation is achieved by directly analysing a raw sequence of timepoints. Frames can then be directly translated to correct drift or a table with drift coordinates can be generated and used in other algorithms (e.g. SRRF).



Drift correction with NanoJ-Core. a) Composite image of two frames from a time-lapse dataset of the same field-of-view. b) Cross-correlation map (CCM) between the two frames shown in a). The vector position of the maximum indicates the linear shift between the two frames. c) Overlay of the two frames after drift correction using NanoJ-Core. d) Vertical and horizontal drift curves obtained using NanoJ-Core from the 100-frame raw data.

NanoJ-Core: Channel Registration


Channel registration is able to calculate and apply an elastic transform to realign different captured wavelengths.



Multi-colour channel registration with NanoJ-Core. a) Composite image of multi-colour beads, prior to (left) and after (right) channel registration using NanoJ-Core. Scale bars: 25 μm , insets: 0.5 μm . b) Vectorial representation of the shift between the two channels (left, displacement vector length 50 times larger for representation purposes), horizontal (middle) and vertical (right) shift maps obtained and applied to the data shown in a). Scale bars: 25 μm .

NanoJ-SRRF: Live-Cell Super-Resolution


Super-Resolution Radial Fluctuations (SRRF) is able to generate a super-resolution reconstruction by analysing fluctuations in the emission of fluorophores captured in a short image sequence.



Live-cell super-resolution microscopy with NanoJ-SRRF. a) Comparison of widefield and SRRF reconstruction from UtrCH-GFP actin labelling. Scale bar: 5 μm . b) Time-course of the inset shown in a), obtained at 33.3 Hz and displayed every 30 s. Scale bar: 1 μm . c) Colour-coded time course. Scale bar: 1 μm .

NanoJ-SQUIRREL: Estimating Image Quality & Resolution


SQUIRREL quantitatively calculates local quality and resolution in super-resolution images. It highlights limitations on the data collected and a reference point to help researchers improve imaging fidelity.



Quality and resolution assessment with NanoJ-SQUIRREL. a) A super-resolution rendering and acquired widefield image of fixed Alexa647 labelled microtubules. b) Left: SQUIRREL error map highlighting discrepancies between the super-resolution and diffraction-limited images in (a). Right: Magnified insets at indicated positions on error map. c) Left: SQUIRREL resolution map of the super-resolution image in (a). Right: Magnified insets for indicated resolution blocks. Whole image scale bars = 5 μm , inset scale bars = 1 μm .

NanoJ-VirusMapper: Structural Mapping and Modelling


VirusMapper features a single-particle analysis (SPA) algorithm, capable of generating structural models of conserved structures imaged by Super-Resolution Microscopy



Quantitative SPA-based modelling. Multicomponent model of the Vaccinia virus by imaging in super-resolution hundreds of fluorescently labelled viruses and modelling their structure through VirusMapper

NanoJ-Fluidics: Liquid Exchange

NanoJ-Fluidics (Pumpy McPumpface) entails the control of a simple inexpensive LEGO-based syringe pump array that automates live-to-fixed imaging and sequential labelling of the sample. Protocols can be run automatically in parallel to the acquisition.



Schematics of the NanoJ-Fluidics system. a) 3D side view of a single syringe pump. b) 2D top view of a syringe pump array (representing 4 pumps out of 128 maximum) and a fluid extraction peristaltic pump, both controlled by an Arduino UNO. c) Example of possible workflows

Author representation

✉ I recognise a gender inequality in the author list (2F/8M). I am committed to achieve a better balanced representation of gender and minorities in our future work.

References

- Almada et al. Automating multimodal microscopy with nanoj-fluidics. *Nature communications*, 10(1):1223, 2019.
- Culley et al. Quantitative mapping and minimization of super-resolution optical imaging artifacts. *Nature methods*, 15(4):263, 2018.
- Gray et al. Virusmapper: open-source nanoscale mapping of viral architecture through super-resolution microscopy. *Scientific reports*, 6:29132, 2016.
- Gustafsson et al. Fast live-cell conventional fluorophore nanoscopy with imagej through super-resolution radial fluctuations. *Nature communications*, 7:12471, 2016.
- Laine et al. Nanoj: a high-performance open-source super-resolution microscopy toolbox. *Journal of Physics D: Applied Physics*, 52(16):163001, 2019.

Funded by

